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Bio-Zertifikat	<input checked="" type="checkbox"/>	
Typ (ISO 22674)	extrahart (4)	
Dichte [g/cm <sup>3</sup> ]	18,8	
WAK [10 <sup>-6</sup> K <sup>-1</sup> ]	25–500 °C 25–600 °C	14,2 14,4
Gießtemperatur [°C]	1270	
Schmelzintervall [°C]	1045–1100	
Elastizitätsmodul [GPa]	100	
Dehngrenze (R <sub>0,2</sub> ) [MPa]	* 420/500	
Zugfestigkeit (R <sub>m</sub> ) [MPa]	570	

Bruchdehnung (A <sub>5</sub> ) [%]	* 10/7
Vickershärte (HV 5)	* 185/215
BEGO-GOLD-Farbcode	gelb 5
Mittlere Korngröße [µm]	20
Entspricht DIN EN ISO 9693	
DIN EN ISO 22674	

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Au	86,0
Pt	11,5
Zn	1,6
Fe, Rh, In	

## Alloy characteristics

Biocertificate	<input checked="" type="checkbox"/>	
Typ (ISO 22674)	extra-hard (4)	
Dichte [g/cm <sup>3</sup> ]	18,8	
CTE [10 <sup>-6</sup> K <sup>-1</sup> ]	25–500 °C 25–600 °C	14,2 14,4
Casting temperature [°C]	1270	
Melting interval [°C]	1045–1100	
Modulus of elasticity [GPa]	100	
Elongation limit (R <sub>0,2</sub> ) [MPa]	* 420/500	
Tensile strength (R <sub>m</sub> ) [MPa]	570	

Ductile yield (A <sub>5</sub> ) [%]	* 10/7
Vickers hardness (HV 5)	* 185/215
BEGO-GOLD colour code	yellow 5
Average grain size [µm]	20
According to DIN EN ISO 9693	
DIN EN ISO 22674	

\* after casting / after firing

## Safety hint

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## Instruction for use

**Modeling/Sprue system:** Minimum wall thicknesses (after finishing): for ceramic veneering 0.4 mm, for acrylic veneering with retention beads 0.3 mm. Connecting bars between pontics as thick and, in particular, as high as possible. Palatinal reinforcement of long bridges with interdental garlands.

**Investing/preheating:** Use phosphate-bonded crown and bridge investment materials (e. g. BellaStar, Bellavest®). Preheating temperature 850 °C.

**Melting/casting:** General: Do not overheat alloy. Use only clean crucibles, one crucible per alloy. Recommendation: to enable an exact identification of each case cast new metal only. In case of re-casting: only re-cast identical alloys. Blast old material. Add at least 50 % of new material. Use ceramic or graphite crucible. Ceramic crucible: sprinkle some Auromelt HF melting powder (REF 52525) on the casting ingots.

**Moment for casting:** Vacuum pressure casting with induction heating (Nautilus®) and centrifugal casting with induction heating (Fornax®). When the last solid component has submerged completely in the melt, continue to heat for 8 to 18 seconds depending on the induction capacity of the casting unit, then trigger casting. Follow operating instructions for Fornax® and Nautilus®. Flame centrifugal casting (Fundor propane/oxygen): As soon as the last solid component has completely submerged in the melting bath, continue to heat for 3 to 5 seconds, then trigger casting. With resistance heating approx. 3 minutes, follow operating instructions of equipment manufacturer.

**After deflasking:** Blast with Korox® 110 at 2–3 bar. Use fine-toothed carbide cutter for finishing.

**Ceramic:** Use veneering ceramics in accordance with DIN EN ISO 9693. Recommendation: Ceramic with reduced firing temperature (e. g. BeCe PRESS [only with powder opaquer], VM 13). Also suitable: Ceramic with firing temperatures up to approx. 960 °C (e. g. IPS Inline, Heraceram, Vintage Halo). Always follow the ceramic manufacturer's instructions. Support frames carefully (use refractory wool, additional wires or make individual firing supports of BegoForm® [REF 52785]). Set the heating rate to a maximum of 55 °C per minute!

**Oxide firing:** Always blast the surface to be veneered (Korox® 110 (REF 46044) at 2 – max. 3 bar) and clean the frame thoroughly (steam clean or boil in aqua dest.). Firing: With vacuum for 5 minutes at 900 °C. Acid treatment of oxide (60 °C, 10–15 minutes). Use a separate bath. Clean thoroughly (brush under running water and then steam clean or boil in aqua dest.).

**Ceramic firing:** Normal cooling is allowed with frames of up to 4 units. Otherwise long-term cooling is recommended.

**Acrylic veneering:** The corresponding instructions of the manufacturers must be followed when processing the veneering systems.

**Final work:** To simplify rubber-polishing, blast-polish with Perblast® micro (lead-free soda glass). After that rubber-polish with BEGO rubber polishers and polish with BEGO precious-metal prepolishing and repolishing compound. Finally clean thoroughly (steam-blast or boil in distilled water).

**Soldering:** Support object in a soldering block of Bellatherm® (REF 51105). Prepare a gap of max. 0.2 mm with parallel walls. Soldering before firing with the flame: PontoStar® G-Solder (REF 61045) and Minoxid (REF 52530) or Fluxsol (REF 52531). Soldering after firing in furnace BEGO-Gold-Solder I (REF 61017) and Minoxid. Allow to cool normally. Acid treatment for residues of Minoxid (60 °C, approx. 1 minute). Clean thoroughly (steam clean or boil in aqua dest.).

**Laser welding:** Filler material: Bio PontoStar® XL wire Ø 0.35 mm (REF 61167).

**Secondary effects:** Such as allergies to contents of the alloy or electrochemically based reactions may very rarely occur.

**Reciprocal actions:** In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

**Reactions:** In case of known incompatibilities and allergies to contents of the alloy.

**Warranty:** Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our products are subject to a constant further development. Therefore alterations in construction and composition are reserved.