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
**Ceramic:** Use ceramics according to ISO 9693 with firing temperatures of up to approx. 980 °C (e. g. BeCe PRESS, VM 13, Omega 900, Duceram Kiss, InLine). Also suitable: Ceramics with reduced firing temperature (e. g. Omega 900). Always follow the ceramic manufacturer's instructions! Always blast the surface to be veneered (pure Korox® 250: 3–4 bar) and clean the frame thoroughly (steam clean or boil in aqua dest.). After cleaning, hold with artery clamps and refrain from touching.

**Oxide firing:** Not necessary. If oxide firing is carried out to check the surface (950–980 °C/5 minutes with vacuum), always blast the oxide again (with pure Korox® 250: 3–4 bar). Clean thoroughly (steam clean or boil in aqua dest.).

**Ceramic firing:** Always apply basic material in two firing operations. The first coating thin (washbrand), the second coating opaque. Wash off frame under running water before application of next ceramic coating. Long-term cooling recommended (cooling phase up to approx. 600 °C). Remove ceramics only mechanically. Hydrofluoric acid (HF) corrodes the metal frame.

**Acrylic veneering:** The corresponding instructions of the manufacturers must be followed when processing the veneering systems.

**Final work:** Blast visible metal surfaces with Korox® 50, then blast-polish external surfaces with Perlblast®. Rubber-polish with BEGO rubber polishers, and after that finish-polish with BEGO cobalt-chrome polishing paste (blue). Clean thoroughly (steam clean or boil in aqua dest.).

**Soldering:** Soldering prior to firing with the flame: Wirobond® solder (REF 52622) and Fluxsol flux (REF 52531). Soldering after firing in the furnace: WGL solder (REF 61079) and Minoxyl flux (REF 52530). Long-term cooling recommended (cooling phase up to approx. 600 °C).

**Laser welding:** Filler material: Wiroweld wire Ø 0.35 mm (REF 50003) or Ø 0.5 mm (REF 50005).

**Prescription only / For professional use only **

**Secondary effects:** Such as allergies to contents of the alloy or electrochemically based reactions may very rarely occur.

**Reciprocal actions:** In case of occlusal or approximal contact of different alloys electrochemically based reactions may very rarely occur.

**Reactions:** In case of known incompatibilities and allergies to contents of the alloy.

**Warranty:** Whether given verbally, in writing or by practical instructions, our recommendations for use are based upon our own experience and trials and can only be considered as standard values. Our products are subject to a constant further development. Therefore alterations in construction and composition are reserved.